



# Autoclavable Nichipet 7000 Multi Channel Digital Pipette

## FEATURES

- Totally Autoclavable. (121°C , 20 minutes)
- UV Resistance.
- 4 sizes covering a wide volume range (0.5-300ul).
- Fatigue free operation with enlarged finger-rest.
- Easy and convenient single hand locking lever.
- Revolvable handle to fix convenient position in pipetting.

### Specifications

Cat. No.	Channel	Volume Range(ul)	Usable Tips or equivalent	Increment
00-NP7-8V	8	0.5-10	BMT-UTR	0.01ul
00-NP7-8S		5-50	BMT-SGR	0.1ul
00-NP7-8L		40-200	BMT-SGR	0.2ul
00-NP7-8K		50-300	BMT-KR	1ul
00-NP7-12V	12	0.5-10	BMT-UTR	0.01ul
00-NP7-12S		5-50	BMT-SGR	0.1ul
00-NP7-12L		40-200	BMT-SGR	0.2ul
00-NP7-12K		50-300	BMT-KR	1ul

### Accuracy / Precision

Cat. No.	Volume	Accuracy(%)	Precision(%)
00-NP7-8V	1ul	±8.0	< 4.0
	5ul	±4.0	< 2.0
00-NP7-12V	10ul	±2.0	< 1.0
	5ul	±3.0	< 1.5
00-NP7-8S	25ul	±2.0	< 0.8
	50ul	±1.0	< 0.5
00-NP7-12S	40ul	±1.4	< 0.5
	100ul	±1.0	< 0.4
00-NP7-8L	200ul	±0.9	< 0.3
	50ul	±1.4	< 0.7
00-NP7-12L	150ul	±1.0	< 0.5
	300ul	±0.6	< 0.2

### Tips(Autoclavable)

Cat. No.	Volume Range(ul)	Color	Usable Model	Q'ty / Box
BMT-UTR	0.1-10	Clear	00-NP7-8V, 12V	960(96pcs. X 10c/s)
BMT-SGR	2-200	Clear	00-NP7-8S, 12S, 8L, 12L	960(96pcs. X 10c/s)
BMT-KR	50-300	Clear	00-NP7-8K, 12K	960(96pcs. X 10c/s)

### Stand

Cat. No.	Storable Number
7000S	2pcs.

### Tray

Cat. No.	Q'ty / Box
88-TY	50

