

Autoclavable Nichipet EX Digital Micro Pipette



FEATURES

- Fatigue free operation with enlarged finger-rest.
- Enhanced UV resistance.
- Light weight and bright color.
- Totally autoclavable.
- 8 sizes cover a wide volume range.(0.1ul-10000ul)
- Easy and convenient single hand lockning lever.
- Filtered nozzles. (1000,5000,& 10000ul)

■ Specifications

Cat. No.	Specifications				
	Volume Range (ul)	Increments (ul)	Measured Volume(ul)	AC (%)	CV (%)
NPX-2	0.1-2	0.002	0.1	*	*
			1	±5.0	<2.5
			2	±3.0	<1.0
NPX-10	0.5-10	0.01	0.5	±5.0	<3.0
			5	±1.0	<1.0
			10	±1.0	<0.5
NPX-20	2-20	0.02	2	±5.0	<3.0
			10	±1.0	<1.0
			20	±1.0	<0.4
NPX-100	10-100	0.1	10	±2.0	<1.0
			50	±1.0	<0.3
			100	±0.8	<0.3
NPX-200	20-200	0.2	20	±1.0	<0.5
			100	±0.8	<0.3
			200	±0.8	<0.2
NPX-1000	100-1000	1.0	100	±1.0	<0.5
			500	±0.8	<0.3
			1000	±0.7	<0.2
NPX-5000	1000-5000	10.0	1000	±1.0	<0.3
			2500	±0.8	<0.3
			5000	±0.6	<0.2
NPX-10ML	1000-10000	10.0	1000	±2.0	<0.4
			5000	±0.8	<0.3
			10000	±0.4	<0.2

* The AC and CV at 0.1ul depend much on the operator's skill.

■ Ordering Information

Usable tips			
Cat. No.	Capacity(ul)	Color	Q'ty/Bag
BMT-UT	0.1-10	Clear	1000
BMT-SS	0.5-10	Clear	1000
BMT-SE	2-200	Clear	1000
BMT-SG	2-200	Clear	1000
BMT-L	100-1000	Blue	1000
BMT-X	1000-5000	Green	200
BMT-Z	1000-10000	Clear	200

(Racked tip is available. Please see page 8.)

Rotary Stand	
MLT-STD	6 Pipettes can be stored.

