

**NEW**

Glass Tip Pipette

# Nichipet ECO

## FEATURES

- Digital micro pipette with GLASS TIP
- GLASS TIP enables totally bias free liquid handling for endocrine disruptors analysis
- 4 sizes covering a wide volume range (20ul to 10000ul)
- Nozzle with double O-rings enables easy tip fitting / ejection and reliable sealing secured
- Nozzle is made of PTFE and piston is made of ceramics

## Specifications

Cat. No.	Volume Range	Incr.	Accuracy	Precision	Usable Tips
00-NPC-200	20-200ul	0.2ul	$\pm < 6.0 - \pm < 3.0$	$< 6.0 - < 3.0$	Glass Tip S
00-NPC-1000	100-1000ul	1.0ul	$\pm < 5.0 - \pm < 3.0$	$< 5.0 - < 3.0$	Glass Tip L
00-NPC-5000	1000-5000ul	10.0ul	$\pm < 5.0 - \pm < 3.0$	$< 5.0 - < 3.0$	Glass Tip X
00-NPC-10ML	1000-10000ul	50.0ul	$\pm < 6.0 - \pm < 3.0$	$< 6.0 - < 3.0$	Glass Tip Z

## Glass Tip

Cat. No.	Volume Range	Package
00-GLT-S	20 to 200ul	10pcs
00-GLT-L	100 to 1000ul	10pcs
00-GLT-X	1000 to 5000ul	10pcs
00-GLT-Z	1000 to 10000ul	10pcs



Further research into " Environmental Hormone " is now strongly requested to clarify the precise cause and effect on the human system and wild living things.

Standard pipette tips made from plastic are not suitable for the research of Environmental Hormone, and the glass is the most reliable and stable material for this research.

