

# Nichipet F&V

Fixed Volume and Triple Volume Micro Pipette

## FEATURES

- 11 different standard volumes for Nichipet F from 5 to 1000 microliter. Special volume in this range is available to order.
- 3 different models, 10/20/50 microliter, 50/100/200 microliter and 200/500/1000 microliter, for Nichipet V
- Smooth movement of the plunger and lightweight. Slender design reduce operator's fatigue.
- Volume change of Nichipet V is simply done by just turning the thumb knob.
- Teflon® seal and O-ring eliminate the need for routine maintenance and lubrication.
- Built-in tip ejector.

### Nichipet V Triple Volume

Three volumes can be selected by turning the thumb knob until the desired volume figures on the thumb correspond to the arrow mark printed on the top of the pipette handle. Selected volume remains securely fixed by clickstop mechanism of the thumb knob.



### Nichipet F Fixed Volume

Most useful 11 fixed volumes are available as standard. Special fixed volume are available to order.



## Ordering Information / Accuracy and Precision

### Nichipet F FIXED VOLUME(Standard Volume) & Nichipet V TRIPLE VOLUME

	Cat. No.	Volume Range(ul)	Accuracy(%)	Precision(%)
V	NPV-S	10/20/50	$\pm 1.2 - \pm 0.7$	$< 0.6 - < 0.3$
	NPV-M	50/100/200	$\pm 0.7 - \pm 0.5$	$< 0.4 - < 0.3$
	NPV-L	200/500/1000	$\pm 0.7 - \pm 0.5$	$< 0.3 - < 0.2$
F	NPF-5	5	$\pm 1.5$	$< 0.5$
	NPF-10	10	$\pm 1.2$	$< 0.4$
	NPF-20	20	$\pm 0.8$	$< 0.3$
	NPF-25	25	$\pm 0.8$	$< 0.3$
	NPF-50	50	$\pm 0.6$	$< 0.3$
	NPF-100	100	$\pm 0.6$	$< 0.3$
	NPF-200	200	$\pm 0.6$	$< 0.2$
	NPF-250	250	$\pm 0.6$	$< 0.3$
	NPF-300	300	$\pm 0.6$	$< 0.3$
	NPF-500	500	$\pm 0.5$	$< 0.2$
	NPF-1000	1000	$\pm 0.5$	$< 0.2$

### Stand

Cat. No.	6 pipettes can be stored.
BM-STD (Acryl resin)	

