



**National Scientific
Supply Company, Inc.**
Plastic Products Manufacturing Division
250 York Place, Claremont, California 91711-4883

Because Quality Matters

(800) 793-1779
Phone: (909) 621-4585
FAX: (909) 621-5653

Certificate of Analysis

DNase/RNase/Endotoxin-Free

This certifies that the following product number(s) conform to all claims stated on the label.

Product Number	Lot Number
RandomSample-NAT-SCI-1	12345
RandomSample-NAT-SCI-2	67890

Products certified RNase and DNase-Free have had randomly selected product assayed for RNase by assessing the integrity of RNA, and for DNase by assessing the integrity of plasmid DNA.

Products certified non-pyrogenic have been tested for bacterial endotoxins. Samples selected at random were tested and validated using the kinetic turbidimetric method with a sensitivity of <.05 EU/mL.

Products certified sterile have been validated for the proper gamma radiation dosage using Method 1 Validation of Radiation Sterilization of the Association for Advancement of Medical Instrumentation (AAMI). This validation procedure is verified to insure all sterilized product receives the adequate doses of gamma radiation to provide a Sterility Assurance Level (SAL) of 10^{-3} . All of our lots are sterilized to the proper dosage levels.

Name: _____

Date: _____

Position: _____